

CEREC Bur Setup

<u>Crown Prep Setup</u>			
	<u>Manufacture</u>	<u>Patterson Order Number</u>	<u>Purpose</u>
Occlusal reduction (depth cuts)	Meisinger	509-9775	Router bur (2mm depth cuts)
Barrel diamond (230 C)	Premier	571-8333	Reduce occlusion to depth cuts
1171	Midwest	493-0053	Break proximal contacts
Flat end taper 2002.8	Premier	572-7193	Axial shoulder (buccal & lingual)
Flat end taper 2001.8	Premier	572-7151	Axial shoulder (mesial & distal)
Flat end taper (fine) 703.8	Premier	572-2707	Smooth margins
Dura-White Stones	Shofu	635-3072	Smooth & Polish prep
<u>Misc Operative</u>			
Inverted cone diamond	Coltene-Whaledent	741-0699	Inlays/Onlays
#4 round	Midwest	492-9550	Caries removal
#6 round	Midwest	492-9568	Caries removal
Razor 330	Axis	120-8545	Metal cutting (amalgam, PFM, etc)
170L	Midwest	492-9840	Proximal boxes
820 EF (skinny)	Meisinger	463-2667	Used to separate interproximals
Finish carbide (egg shape)	Midwest	493-7728	Finish
Finish flame (T2)	Midwest	493-7835	Finish
Finish football diamond	NTI	545-3329	adjust emax